EA020 ELECTROSURGICAL GENERATOR USER MANUAL

Important Warranty Information:

Limited Warranty. Applied Medical warrants to the original purchaser for a period of one (1) year from the date of shipment that the generator ("product") will: (i) conform to the then-current written product specifications which may be amended from time to time; and (ii) be free from defects in material and workmanship under normal use and preventative maintenance.

Troubleshooting. To return the product, the purchaser must first contact Applied Medical's customer service department for troubleshooting by Applied Medical. After such contact, Applied Medical may provide authorization to return ship the product to Applied Medical.

Product Return. Upon receipt of authorization, the product may be returned freight prepaid or by Applied Medical field representatives. Any unauthorized or freight collect returns will be refused and/or returned to the purchaser at the purchaser's expense. Used product returned to Applied Medical must be cleaned and packaged for safe handling. Applied Medical's field representatives are prohibited from picking up product that does not meet these criteria. In the event that the product is repaired or replaced by Applied Medical, the warranty period will not be extended.

Warranty Claims. Applied Medical's sole obligation under the limited warranty set forth above, is for Applied Medical to, at its option: (i) repair or replace the product or any part thereof at its sole expense; or (ii) issue a credit for the net invoiced purchase price, provided that the product is returned to Applied Medical during the applicable time period of the limited warranty and, is determined by Applied Medical to be defective.

Exclusions. Applied Medical's obligations under this limited warranty are subject to purchaser's compliance with all instructions and requirements regarding the use of the product. This limited warranty does not cover any product that has been (1) improperly used or stored, or subjected to physical abuse such as dropping, negligence, recklessness, willful misconduct, or any failure to follow the procedures set forth in the product manual or as demonstrated by Applied Medical representatives; (2) repaired or modified by any person or entity other than Applied Medical or not authorized in writing by Applied Medical; or (3) damaged resulting from the use of non-Applied Medical devices, components or parts. This limited warranty shall be personal to the original purchaser. Any sale, rental or other transfer or use of the product covered by this warranty to or by a user other than the original user shall cause this warranty to immediately terminate. This limited warranty does not apply to any accessories or devices used with the product. Applied Medical reserves the right to make changes to the product at any time without creating an obligation to make the same or similar changes to products previously built or sold by Applied Medical.

Force Majeure Events. Applied Medical will not be liable for damages or product failures or failure to perform or delay due to acts of God, strikes, fires, explosions, flood, riot, unavoidable accidents, power disturbances, acts of government, or similar occurrences outside of the control of Applied Medical.

DISCLAIMER. EXCEPT AS PROVIDED EXPRESSLY ABOVE OR TO THE EXTENT PROHIBITED BY APPLICABLE LAW, APPLIED MEDICAL EXCLUDES ALL WARRANTIES, EXPRESS OR IMPLIED, BY OPERATION OF LAW OR OTHERWISE, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. APPLIED MEDICAL'S LIABILITY WITH RESPECT TO THE LIMITED WARRANTY SHALL BE LIMITED TO THE PURCHASE PRICE FOR THE PRODUCTS SOLD TO PURCHASER. IN NO EVENT SHALL APPLIED MEDICAL BE LIABLE FOR ANY SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, OR CONSEQUENTIAL LOSS, DAMAGE, OR EXPENSE, HOWEVER CAUSED, ON ANY THEORY OF LIABILITY (INCLUDING LOST PROFITS, PROCUREMENT OF SUBSTITUTE GOODS, AND BUSINESS INTERRUPTION), AND WHETHER OR NOT APPLIED MEDICAL HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES. APPLIED MEDICAL NEITHER ASSUMES, NOR AUTHORIZES, ANY THIRD PARTY TO ASSUME FOR IT, ANY LIABILITY OR RESPONSIBILITY IN CONNECTION WITH THE PRODUCT.

Software License Agreement:

LICENSE. The software embedded within the Applied Medical generator ("System") and the accompanying documentation (collectively, "Software") are provided under license and are not sold to purchaser. Applied Medical Resources Corporation ("Applied Medical") grants to purchaser a limited, non-exclusive, non-sublicensable, non-transferable, fully-paid, revocable license to (a) use the Software solely as incorporated in the System in machine-executable object code form on a single central processing unit embedded in the System, for internal business purposes, and solely in connection with the operation of the System as described in the documentation, and (b) use the documentation solely for the purpose of using the Software in compliance with this license.

Purchaser shall not (and shall not cause or allow a third party to) (i) use, copy, translate, modify, create derivative works of (including without limitation, translations, transformations, adaptations, or altered versions), or transfer, (ii) merge the Software with any other software or product not supplied by Applied Medical, (iii) sublicense, lease, rent, loan, sell, or otherwise transfer, (iv) reverse engineer, decompile, disassemble, attempt to derive the source code for, or otherwise manipulate, or (v) distribute, disclose, permit to be disclosed or publicly display or perform, in any forum or by any means, the Software, in whole or in part, or any copy thereof. Manipulation of the Software is permitted only to the extent the Software is required to be modified by applicable law, but only upon prior 30 days written notice to Applied Medical.

Purchaser shall not copy the Software under any circumstances. Purchaser shall not alter, obscure or remove any copyright, trademark, proprietary rights, disclaimer, or warning notice included on or embedded in any part of the Software (including those of third parties).

OWNERSHIP. Except for the limited rights expressly granted herein, Applied Medical and its affiliates own the entire right, title and interest to the Software and any modifications, improvements, and derivations thereof. The Software is licensed, not sold, to purchaser. There is no implied license, right or interest granted in any copyright, patent, trade secret, trademark, service mark, or other intellectual property right. Nothing in the limited license requires Applied Medical to produce or furnish technical data to purchaser.

WARRANTY. Any warranty for the System is included with the system documentation.

DISCLAIMER. EXCEPT AS PROVIDED IN THE SYSTEM WARRANTY OR PROHIBITED BY APPLICABLE LAW, APPLIED MEDICAL EXCLUDES ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, BY OPERATION OF LAW OR OTHERWISE, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT SHALL APPLIED MEDICAL BE LIABLE FOR ANY SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, OR CONSEQUENTIAL LOSS, DAMAGE, OR EXPENSE, HOWEVER CAUSED, ON ANY THEORY OF LIABILITY (INCLUDING LOST PROFITS, PROCUREMENT OF SUBSTITUTE GOODS, AND BUSINESS INTERRUPTION), AND WHETHER OR NOT APPLIED MEDICAL HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES. APPLIED MEDICAL NEITHER ASSUMES, NOR AUTHORIZES, ANY THIRD PARTY TO ASSUME FOR IT, ANY LIABILITY OR RESPONSIBILITY IN CONNECTION WITH THE SOFTWARE.

TERM. This license will begin on the date the Software is delivered to purchaser. Notwithstanding the foregoing, this license shall terminate immediately upon written notice to purchaser by Applied Medical if purchaser breaches any term or condition of this license. Purchaser agrees upon notice of termination to promptly discontinue all use of the Software and destroy all copies thereof (whether in tangible form or as installed on purchaser's equipment).

EXPORT LAWS AND REGULATIONS. The Software and related technology may be subject to U.S. export control laws and export or import regulations in other countries. Purchaser shall strictly comply with all such laws and regulations and has the responsibility to obtain such licenses to export, re-export or import as may be required.

U.S. GOVERNMENT PURCHASERS. The Software is a "commercial item" as that term is defined at 48 C.F.R. 2.101, developed exclusively at private expense, consisting of "commercial computer software" and "commercial computer software documentation" as such terms are used in applicable U.S. acquisition regulations, all U.S. government end users acquire the Software with only those rights set forth therein.

MISCELLANEOUS. This software license sets forth the entire agreement and understanding of the parties relating to the Software and supersedes all prior agreements, understandings, and discussions between them relating to the Software. No modification of or amendment to this software license, nor any waiver of any rights under this software license, shall be effective or binding on the parties unless made in writing and signed by an authorized representative of each of the parties. No failure on the part of a party in exercising any power or right hereunder shall operate as a waiver thereof; nor shall any single or partial exercise of any power or right preclude any other or further exercise thereof or the exercise of any other power or right.

Whenever possible, each provision of this software license shall be interpreted in such a manner as to be effective and valid under applicable law. However, in the event that any provision becomes prohibited or invalid under applicable law, or is otherwise held unenforceable, then such provision shall be modified to reflect the parties' intent, consistent with applicable law invalidating the remainder of such provision or the remaining provisions of the software license.

Third-Party Software Acknowledgement:

Portions of the graphics display use emWin, which is licensed under a Single Product Source Code License (SP-SCL) with SEGGER.

Voyant Electrosurgical Generator User Manual

Manufactured By
APPLIED MEDICAL RESOURCES CORP.
22872 Avenida Empresa
Rancho Santa Margarita, California 92688 USA
TEL: 949-713-8300
FAX: 949-713-8400
www.appliedmedical.com

European Authorized Representative
Applied Medical Europe BV
Wiekenweg 21
3815 KL Amersfoort
The Netherlands

Australian Sponsor
Applied Medical Australia
1/32 Windorah Street
Stafford QLD 4053 Australia

Patents: patents.appliedmedical.com

Applied Medical and the Applied Medical logo design are U.S. registered trademarks of Applied Medical Resources Corporation.

©2018 Applied Medical Resources Corporation. All rights reserved.

Made in the USA

Prices, specifications, and model availability are subject to change without notice.

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Voyant Electrosurgical Generator User Manual

EA020 Electrosurgical Generator User Manual	5
Manual del usuario del generador electroquirúrgico EA020	33
EA020 manuel d'utilisation du générateur électrochirurgical	61
Manuale utente del Generatore elettrochirurgico EA020	89
Elektrochirurgie-Generator EA020 – Benutzerhandbuch	117
Gebruikershandleiding van de EA020 elektrochirurgische generator	145
EA020 elektrokirurgisk generator - användarhandbok	173
Brugsvejledning til EA020 elektrokirurgisk generator	201
Sähkökirurgisen EA020-generaattorin käyttöopas	229
EA020 bruksanvisningen for elektrokirurgisk generator	257
EA020 Manual do Utilizador do Gerador Eletrocirúrgico	285

EA020 ELECTROSURGICAL GENERATOR USER MANUAL

Chapter 1: Introduction	8
Voyant Electrosurgical Generator Overview	
Indications for Use	
Product Contents	
Controls & Receptacles	
Chapter 2: Patient & Operating Room Safety	10
Terminology Used	
General Warnings & Safety Instructions	
Fire/Explosion Warnings & Safety Instructions	
Electromagnetic Compatibility (EMC) Warnings & Safety Instructions	
Chapter 3: Generator Setup & Operation	12
Setting Up the Voyant Electrosurgical Generator	
Powering On the Voyant Electrosurgical Generator	
Connecting Voyant Devices	
Sealing Tissue with Voyant Devices	
Alarm Conditions	
After Surgery	
Chapter 4: Generator Options & Info	16
Navigating the Main Menu	
Adjusting the Generator Tone Volume	
Adjusting the Display Brightness	
Accessing Generator Information	
Changing the Display Language	
Chapter 5: Troubleshooting & System Alarms	17
Troubleshooting	
System Alarms	
Chapter 6: Cleaning, Maintenance, Service & Disposal	20
Cleaning	
Maintenance	
Service	
Disposal	
Chapter 7: Technical Specifications	24
Conformance to Standards	
Generator Specifications	
Power Cable Specifications	
Protection against Electric Shock	
Equipotential Stud	
RF Energy Output Data	
Electromagnetic Compatibility	
Explanation of Symbols	313

Voyant Electrosurgical Generator Overview

The Voyant electrosurgical generator (model EA020) outputs radiofrequency (RF) energy to Voyant devices for the sealing of vessels and tissue bundles. The Voyant generator is a reusable power supply that operates outside of the sterile field.

Indications for Use

The Voyant electrosurgical generator is indicated for use with Voyant devices in open and laparoscopic surgical procedures where the electrosurgical ligation of vessels or tissue bundles is desired.

Note: For indications specific to each Voyant device used with the Voyant generator, refer to each device's Instructions for Use (IFU).

Product Contents


The shipping box contains the Voyant generator, power cable, verification adapter, and user manual. The EX0xx series power cables, EX060 equipotential cable, EX240 verification adapter, and Voyant cart are also available separately.

Controls & Receptacles

Voyant Electrosurgical Generator – Front Panel

Voyant Electrosurgical Generator – Back Panel

- 1. Stand-By Button
- 2. Soft Keys

- 3. Display
- 4. Device Port

- 5. Power Entry Module
- 6. Equipotential Stud

Terminology Used

- **WARNING:** Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury to the patient/user.
- **PRECAUTION:** Indicates a potentially hazardous situation which, if not avoided, could result in moderate or minor injury to the patient/user or cause undesirable device effects.

Note: Indicates information that is important, but not directly related to patient/user safety.

General Warnings & Safety Instructions

WARNING – To reduce the risk of death or serious injury, observe the following warnings:

- **WARNING:** The generator and associated devices are for use by qualified medical personnel only. Do not use the system unless properly trained on the use of the system, electrosurgery, and the specific procedure to be performed.
- **WARNING:** Read all instructions, warnings, precautions, and notes provided in this manual prior to operating the generator. Refer to the Instructions for Use (IFU) provided with the Voyant devices for device-specific instructions, warnings, precautions, and notes.
- **WARNING:** Carefully inspect the generator, power cable, and devices for any signs of damage prior to use. Do not use any damaged equipment.
- **WARNING:** Do not connect any non-Voyant devices or power cables other than those specified for the generator. See Chapter 7, Technical Specifications, for power cable specifications.
- **WARNING:** Use electrosurgical devices with caution in the presence of cardiac pacemakers or other active implants (e.g. implantable cardioverter defibrillators). Electrosurgical devices can interfere with the function of pacemakers and other active implants. To ensure patient safety, obtain qualified advice from the pacemaker or active implant manufacturer, or hospital Cardiology Department, prior to using electrosurgical devices on patients with cardiac pacemakers or other active implants.
- **WARNING:** The patient should not come into contact with metal objects which are earthed or which have an appreciable capacitance for earth (e.g. operating table supports or I.V. poles). The use of antistatic sheeting is recommended.
- **WARNING:** Do not activate more than one electrosurgical device at a time as unintended tissue effects and/or product damage may result.
- **WARNING:** When electrosurgical devices and physiological monitoring equipment are used simultaneously on the same patient, any monitoring electrodes should be placed as far as possible from the surgical site. Needle monitoring electrodes are not recommended. Monitoring systems incorporating high-frequency current limiting devices are recommended when electrosurgical devices are used.
- **WARNING:** System failure could result in an unintended increase or decrease in output power, a change in the output waveform, or the cessation of electrosurgical output. A backup system/method should be available for use.

PRECAUTION – To reduce the risk of moderate to minor injury and unintended surgical effects, observe the following cautions:

- **PRECAUTION:** High-frequency surgical equipment can produce a physiological effect, such as neuromuscular stimulation.
- **PRECAUTION:** Studies have shown that the smoke generated by electrosurgical devices can be potentially harmful to the patient/user. Use a smoke evacuation device when necessary.
- **PRECAUTION:** Only use instruments that can withstand the maximum output (rms) voltage for each output mode (refer to specific device IFU for rated accessory voltage). Using a voltage range that is greater than the maximum output voltage may result in patient or operator injury, or may damage the instrument.

Note: The mains plug or appliance coupler may be used as the disconnect device. Position the generator such that either end of the power cable is readily available.

Fire/Explosion Warnings & Safety Instructions

WARNING – To reduce the risk of fire, explosion, death or serious injury, observe the following warnings:

- **WARNING:** Allow flammable agents used for cleaning or disinfecting or as solvents for adhesives to evaporate before electro-surgical devices are used. Use non-flammable agents for cleaning and disinfection whenever possible.
- **WARNING:** There is a risk of flammable solutions pooling under the patient or in the patient's body depressions/ cavities. Any pooled fluid should be removed before electro-surgical devices are used.
- **WARNING:** Do not use electro-surgical devices in the presence of flammable anesthetics. Verify that all anesthesia circuit connections are leak free before and during the use of electro-surgical devices.
- **WARNING:** Avoid the use of oxidizing gases such as nitrous oxide (N₂O) and oxygen (O₂) near electro-surgical devices.
- **WARNING:** Prevent flammable or oxidizing gases from pooling under surgical drapes or near the surgical site.
- **WARNING:** Avoid the accumulation of flammable endogenous gases.
- **WARNING:** When saturated with oxygen, some materials (such as gauze, cotton or surgical drapes) may be ignited by the sparks and heat that can be produced by the use of electro-surgical devices. Do not place electro-surgical devices near or in contact with flammable materials. Keep gauze wet if it is near electro-surgical devices.

Electromagnetic Compatibility (EMC) Warnings & Safety Instructions

WARNING – To reduce the risk of death or serious injury, observe the following warnings:

- **WARNING:** Interference produced by the operation of high-frequency surgical equipment, such as electro-surgical generators/devices, may adversely influence the operation of other electrical equipment. Ensure that the generator is placed an adequate distance from other electronic equipment. Do not cross or bundle electronic device cords.
- **WARNING:** Use of accessories and cables other than those specified may result in increased electromagnetic emissions or decreased electromagnetic immunity.
- **WARNING:** No customer modification of this equipment is allowed. Modification of this equipment could have a negative impact on electrical safety and electromagnetic emissions.
- **WARNING:** Avoid use of the generator adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the system should be observed to verify normal operation in the configuration in which it will be used.
- **WARNING:** Electrostatic discharge (ESD) on hand device port pins could cause the generator screen to become blank. Restarting the generator will restore system into normal condition.

PRECAUTION – To reduce the risk of moderate to minor injury and unintended surgical effects, observe the following cautions:

- **PRECAUTION:** The generator requires special precautions regarding EMC and must be installed and used according to the EMC information provided in this manual.
- **PRECAUTION:** Portable and mobile RF communications equipment can affect medical electrical equipment.


Setting Up the Voyant Electrosurgical Generator

- **WARNING: To avoid the risk of electric shock, the generator must only be connected to a supply mains with protective earth. Do not use power plug adaptors, power strips or extension cords.**
 - PRECAUTION: Allow the generator surface temperature to equalize to ambient temperature before turning on the generator.
 - PRECAUTION: Do not stack equipment on top of the generator or place the generator on top of other equipment as the generator may become unstable or be unable to cool properly.
 - PRECAUTION: To avoid overheating of the generator, keep the air vents clear.
 - PRECAUTION: Confirm that the electrical installation of the room complies with local codes and regulatory requirements.
 - PRECAUTION: To avoid the unintentional disconnection of the generator's power cable, place the generator such that accidental contact or access to the power cable is limited.
1. Place the generator on the Voyant cart or a stable surface. To secure the generator on the cart, screw the rubber feet through the holes on top of the cart.
 2. Provide adequate space for ventilation around the back of the generator.
 3. Plug the supplied power cable into the appropriate receptacle on the back of the generator. Plug the other end of the cable into a properly grounded receptacle. Once connected, the stand-by button on the front panel of the generator will illuminate orange.
 4. If required by local codes, connect the equipotential stud on the back of the generator to the hospital equalization connector using an equipotential cable.

Powering On the Voyant Electrosurgical Generator


1. Press the stand-by button to turn on the generator. Once the stand-by button has been pressed, the following will occur:
 - The stand-by button will turn green, and the navigation buttons and Applied Medical logo will illuminate.
 - The generator will run through an initial power on self test (POST), and the screen will display the Voyant logo and a status bar to indicate activity (Figure 1).
 - The generator will emit a short start-up tone upon completion of the POST.

Figure 1:

2. If no device is connected at this point, the generator screen will instruct the user to connect a device (Figure 2).


Figure 2:

Connecting Voyant Devices


- **WARNING: Do not connect a wet device to the generator as electrical shock may occur.**
 - **WARNING: Confirm that the device is securely connected to the generator prior to use. Ensure that no metal is exposed at any connection points.**
 - **WARNING: Keep device cords away from the cables of other electrical equipment. Do not wrap device cords around metal objects. Position device cords in such a way that patient contact is avoided.**
 - **PRECAUTION: Do not insert fingers or objects other than the correct connector into the receptacle of the generator.**
Note: When using bipolar Voyant devices with the generator, return pads are not necessary.
1. Connect a Voyant device to the device port on the front of the generator.
 - When a Voyant device is connected, the following will occur:
 - The generator will emit a tone upon initial device connection and the generator screen will momentarily display “Loading Device” (Figure 3).

Figure 3:

- Once the device is ready for use, the rectangle on the generator screen will turn green (Figure 4).

Figure 4:

Sealing Tissue with Voyant Devices

- **WARNING:** Prior to use, refer to the device IFU for the indicated vessel sealing size.
 - **WARNING:** If the tissue effect is not as expected, do not continue to use the generator and/or device until the concern is resolved.
 - PRECAUTION: Confirm proper generator settings prior to activating the device.
 - PRECAUTION: Ensure that the generator activation tone is audible during the procedure.
Note: The Voyant generator has no user-adjustable Fuse level settings.
1. Activate the Fuse mode by pressing and holding the Fuse button on the device handle.
 - While energy is being delivered, an activation tone will sound and the rectangle on the generator screen will turn purple (Figure 5).

Figure 5:

2. When the seal cycle is complete, an end tone will sound and the generator will automatically stop delivering energy. Release the Fuse button after the seal cycle end tone is heard.
 - The rectangle on the generator screen will return to green after the seal cycle completes.

Alarm Conditions

If an alarm condition occurs, the generator will:

- emit a short, pulsed tone,
- display an alarm message, and
- interrupt energy delivery until the condition is corrected.

See Chapter 5, *Troubleshooting & System Alarms*, for more information.


After Surgery

- PRECAUTION: When disconnecting the power cable, always grasp the cable by the plug. Do not pull on the cable itself.
1. Press the stand-by button on the front panel of the generator. The stand-by button will turn orange.
 2. Unplug the power cable from the wall outlet.
 3. Carefully detach the device from the generator by depressing the latch on the top of the device key and gently pulling the device key from the receptacle.
 - Discard devices after use according to hospital policy for the disposal of hazardous medical waste.

Navigating the Main Menu

1. To access the main menu, press the soft key that appears underneath the word “Menu” on the main device screen.
2. To select an option, use the soft keys to highlight and select the desired option (Figure 6).

Figure 6:

3. To exit the main menu, press the soft key that appears underneath the word “Exit” on the main menu screen.

Adjusting the Generator Tone Volume

1. Press the soft key that appears underneath the word “Volume” on the main device screen, or access the main menu and select “Volume” from the menu options.
2. Use the soft keys to adjust the volume of the generator tones.

Note:

- The generator tones cannot be turned off completely.
 - The volume of the generator alarms cannot be adjusted separately from the generator tones.
3. To save your changes and exit the volume screen, press the soft key that appears underneath the word “Exit”.

Adjusting the Display Brightness

1. Access the main menu and select “Brightness” from the menu options.
2. Use the soft keys to adjust the brightness of the display.
3. To save your changes and exit the brightness screen, press the soft key that appears underneath the word “Exit”.

Accessing Generator Information

1. Access the main menu and select “Generator Info” from the menu options.
2. To exit the generator info screen, press the soft key that appears underneath the word “Exit”.

Changing the Display Language

1. Access the main menu and select “Language” from the menu options.
2. Use the soft keys to highlight and select a desired language.
3. After selecting a language, the generator will automatically display all text in the selected language.

Note: The selected language will display each time the generator is turned on.

4. To exit the language screen, press the soft key that appears underneath the word “Exit”.

Troubleshooting

- **WARNING: If there is a reoccurring problem, do not use the system. Contact your sales representative or Applied Medical Device Support for service.**

Situation	Suggestion
The generator screen is blank	<ul style="list-style-type: none"> ● Ensure that the generator was turned on via the stand-by button on the front panel ● Check the power cable and wall outlet ● Restart the generator ● Contact your sales representative or Applied Medical Device Support for service, if necessary
The stand-by button is green but the generator screen is blank	<ul style="list-style-type: none"> ● Remove the system from service and contact your sales representative or Applied Medical Device Support
Generator displays the "Connect A Device" message when a device is attached	<ul style="list-style-type: none"> ● Verify that the device key is fully inserted into the generator port ● Check the device key and the generator port for damage ● Replace the device, if necessary
The generator displays an orange alarm message	<ul style="list-style-type: none"> ● See the <i>System Alarms</i> section of this chapter for more information

System Alarms

When an alarm condition occurs, the generator will emit a pulsed tone and the screen will display one of the following messages with no delay. Energy delivery will be interrupted until the condition is corrected.

Low Priority Alarms

Low priority alarms indicate that operator awareness is required. When low priority alarm conditions occur, the screen shows an orange alarm message with a low priority alarm symbol, and a two-pulsed tone will sound.

Medium Priority Alarms

The Voyant generator does not have any medium priority alarms.


High Priority Alarms

The Voyant generator does not have any high priority alarms.

Note: Position the generator so that alarms shown on the generator display are visible to the surgeon or nursing staff from the surgical field.

Note: The Power On Self Test (POST) verifies that the alarm system and audio and visual output are operational. Additional testing to confirm alarm system operation is not required.

Figure 7:

**Meaning**

This may mean that:

- The device jaws are covered in excessive eschar
- An insufficient amount of tissue is between the jaws
- The device jaws are in contact with metal (e.g. staples, clips, etc.)
- The device jaws are submerged in pooled fluids

Solution

Correct the condition (i.e. clean the jaws, grasp thicker tissue, check for metal objects, and/or remove excess fluid) and then regrasp the tissue and reactivate the device. If the message persists, replace the device.

Figure 8:

**Meaning**

The connected device is unrecognized by the generator.

Solution

Verify that the device key is fully inserted into the generator port. If the message persists, replace the device.

Figure 9:


**Meaning**

- The seal cycle is interrupted before the cycle is complete, or
- The maximum seal cycle time is reached

Solution

Reactivate the seal cycle without repositioning the device.

Figure 10:


**Meaning**

The Voyant device's fuse activation button is released by the user before the seal cycle is complete.

Solution

Reactivate the seal cycle without repositioning the device.

Figure 11:

Meaning

A device button may be stuck.

Solution

Check for a stuck device button. Replace the device if necessary.

Figure 12:


Meaning

Internal diagnostics have identified an internal error that the user cannot correct.

Solution

Turn the power off, then back on again. If the message persists, note the error code provided and contact your sales representative or Applied Medical Device Support for service.

Figure 13:

Meaning

The generator is overheating.

Solution

Check for blocked air vents. Wait 5 minutes before turning the power off. Turn the power back on again. If the message persists, note the error code provided and contact your sales representative or Applied Medical Device Support for service.

Figure 14:

Meaning

Internal diagnostics have identified an internal error during Calibration Verification testing.

Solution

Turn the power off, then back on again. If the generator successfully passes the Power On Self Test (POST), access Service Mode and run Calibration Verification again. If the error message persists, contact your sales representative or Applied Medical Device Support for service.

Cleaning

- PRECAUTION: Do not use caustic, corrosive, or abrasive cleaning solutions or materials.
 - PRECAUTION: While cleaning the generator, do not allow fluids to enter the generator or connectors.
1. Turn off the generator and unplug the power cable from the wall.
 2. Thoroughly wipe all surfaces of the generator and power cable following the cleaning procedures approved by your institution. Mild cleaning solutions or disinfectants, such as 70% Isopropyl alcohol should be used.

Maintenance

A periodic safety check is recommended to verify functionality of the system. The periodic safety check includes:

- Inspect the system and power cord for physical damage
- Power On Self Test (POST)
- Verify the audio
- Verify the calibration using the “Calibration Verification” test series described below
- Electrical safety test in accordance with IEC 60601-1

Note: Verification of generator calibration is recommended to be performed every twenty-four (24) months using the EX240 verification adapter. The recommended method for this testing is the “Calibration Verification” mode accessible through “Service Mode” screen. Please see the “Calibration Verification” section below for further details.

Required Equipment:


- EX240 Verification Adapter (provided) (Figure 15)

Figure 15:

Note: When the verification adapter is connected to the generator, the “Service Mode” screen will be displayed (Figure 16). To exit “Service Mode”, disconnect the verification adapter.

Figure 16:


Calibration Verification

The EA020 Voyant Generator is capable of self-verification using the “Calibration Verification” mode found on the Service Mode screen.


1. Connect the verification adapter to the device port.
2. Use the soft keys to highlight and select “Calibration Verification”.
3. Press the soft key located below “Run” to start testing (Figure 17).

Figure 17:


4. The generator will run ten (10) tests to verify its calibration (Figure 18).

Figure 18:

5. Once the generator has completed the tests, detailed results will be shown on the generator display. Use the soft keys to navigate between the results for each of the ten (10) tests (Figure 19). Record the results of each test in Table 1.

Figure 19:

6. Upon exiting the detailed results, a message will appear on the display indicating the outcome of Calibration Verification testing (Figure 20).

Note: If calibration verification testing fails, the generator will emit an alarm tone and display an error message. For instructions on how to resolve this condition, refer to Chapter 5, *Troubleshooting & System Alarms*.

Figure 20:


Table 1: Calibration Verification Test Data

Test	Measurement	Range	Result	Pass/Fail
1	Voltage (V)	174.51 – 213.29		
2	Voltage (V)	72.90 – 80.57		
3	Voltage (V)	54.68 – 60.43		
4	Voltage (V)	54.81 – 60.58		
	Current (mA)	2879.07 – 3182.13		
	Power (VA)	166.12 – 183.60		
	Phase	(-6.37) – (-2.37)		
5	Voltage (V)	18.26 – 20.19		
	Current (mA)	946.58 – 1046.22		
	Power (VA)	18.20 – 20.11		
	Phase	(-6.56) – (-2.56)		
6	Voltage (V)	36.88 – 40.76		
	Current (mA)	1931.55 – 2134.87		
	Power (VA)	74.99 – 82.88		
	Phase	(-6.43) – (-2.43)		
7	Voltage (V)	36.70 – 40.56		
	Current (mA)	1820.66 – 2012.31		
	Power (VA)	70.33 – 77.74		
	Phase	15.11 – 19.11		
8	Voltage (V)	19.29 – 21.32		
	Current (mA)	959.19 – 1060.16		
	Power (VA)	19.48 – 21.53		
	Phase	15.09 – 19.09		
9	Voltage (V)	37.93 – 41.92		
	Current (mA)	1884.44 – 2082.80		
	Power (VA)	75.23 – 83.15		
	Phase	15.08 – 19.08		
10	Impedance (Ohms)	10.00 – 30.00		

Guided RF Output Verification

"Calibration Verification" mode is the recommended method of performing RF output verification testing of the Voyant Electrosurgical Generator. The generator also performs an internal verification test at POST. If hospital policy requires that the generator is also tested with equipment external to the generator, please contact Applied Medical Device Support.

Service

- **WARNING: The Voyant generator has no user-serviceable parts. Do not attempt to remove the cover or perform service on the generator.**

For service or repair, please contact your sales representative or Applied Medical Device Support.

Applied Medical Device Support Contact Information

Australia
Tel: 1.800.666.272
Email: customerservice.australia@appliedmedical.com

Belgium
Tel: 0800 950 14
Email: customerrelations-be@appliedmedical.com

Finland
Tel: 0800 949 16
Email: customerrelations-fi@appliedmedical.com

France
Tel: 0800 11 12 06
Email: customerrelations-fr@appliedmedical.com

Italy
Tel: 800 210 155
Email: customerrelations-it@appliedmedical.com

Luxembourg
Tel: 0800 225 19
Email: customerrelations-be@appliedmedical.com

Norway
Tel: 800 690 86
Email: customerrelations-no@appliedmedical.com

Sweden
Tel: 0200 120 403
Email: customerrelations-se@appliedmedical.com

The Netherlands
Tel: 0800 0200 144
Email: customerrelations-ne@appliedmedical.com

United States
Tel: 800.417.3218
Email: Medsupport@appliedmedical.com

Austria
Tel: 0800 295 840
Email: customerrelations-at@appliedmedical.com

Denmark
Tel: 80 25 14 66
Email: customerrelations-dk@appliedmedical.com

Germany
Tel: 0800 0347 333
Email: customerrelations-de@appliedmedical.com

Ireland
Tel: 1800 94 88 59
Email: customerrelations-ie@appliedmedical.com

New Zealand
Tel: 0800.644.344
Email: customerservice.nz@appliedmedical.com

Portugal
Tel: 800 780 441
Email: customerrelations-pt@appliedmedical.com

Spain
Tel: 900 993 143
Email: customerrelations-es@appliedmedical.com

Switzerland
Tel: 800 56 3807
Email: customerrelations-ch@appliedmedical.com

United Kingdom
Tel: 0800 8766 882
Email: customerrelations-uk@appliedmedical.com

For more contact information, visit www.appliedmedical.com

Disposal

Do not dispose of the Voyant generator or accessories in the trash with other waste. Local, regional, and/or facility disposal processes for electronic waste should be followed. Contact your sales representative or Applied Medical Device Support to discuss waste disposal solutions and processes.

Conformance to Standards

EN (IEC) 60601-1 (with Canadian and US National Deviations)

Medical electrical equipment – Part 1: General requirements for basic safety and essential performance

EN (IEC) 60601-2-2

Medical electrical equipment – Part 2-2: Particular requirements for the basic safety and essential performance of high frequency surgical equipment and high frequency surgical accessories

EN (IEC) 60601-1-2

Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance – Collateral standard: Electromagnetic compatibility – Requirements and tests

EN (IEC) 60601-1-8

Medical electrical equipment – Part 1-8: General requirements for basic safety and essential performance – Collateral standard: General requirements, tests, and guidance for alarm systems in medical electrical equipment and medical electrical systems

Generator Specifications

General	
Generator Weight	14.5 lb (6.6 kg)
Generator Depth	13.8 in (35.1 cm)
Generator Width	12.0 in (30.5 cm)
Generator Height	4.4 in (11.3 cm)
Cooling	Fan-cooled

Operating Parameters	
Ambient Temperature Range	10° to 40°C (50° to 104°F)
Relative Humidity	20% to 75% non-condensing
Atmospheric Pressure	700 to 1060 mBar

Transport and Storage	
Ambient Temperature Range	-20° to 55°C (-4° to 131°F)
Relative Humidity	10% to 95% non-condensing
Atmospheric Pressure	500 to 1060 mBar

Activation Tone	
Generator Tone Volume	48 dB to 66 dB
Frequency and Duration	Fuse: 500 Hz (interrupted tone) followed by 635 Hz end-tone (750 ms)

Alarm Tone	
Alarm Volume	60 dB to 68 dB
Frequency and Duration	Low Priority Alarms: Two pulse tone; 464 Hz (196 ms), 368 Hz (215 ms)

Output Waveform Characteristics	
Waveform	The RF output is a variable amplitude sinusoidal waveform.
Crest Factor	A constant crest factor of 1.4 nominal for all RF outputs.
Power Rating	29 VA to 375 VA, ± 18.75 VA
Normal Voltage Rating	$15 V_{rms}$ to $150 V_{rms}$, $\pm 7.5 V_{rms}$
High Voltage Rating	$150 V_{rms}$ to $300 V_{rms}$, $\pm 30 V_{rms}$
Normal Voltage Mode Current Rating	$0.5 A_{rms}$ to $8 A_{rms}$, $\pm 0.25 A_{rms}$
High Voltage Current Rating	$0.5 A_{rms}$ to $4 A_{rms}$, $\pm 0.25 A_{rms}$
Mode of Operation	Non-continuous; duty cycle: 10 seconds on, 30 seconds off

Note: Voltage limits are determined by the connected Voyant device. Refer to each Voyant device's Instructions for Use (IFU).

Electrical	
Generator Power Source (Operating Range)	100-240 V, 50-60 Hz, 7-3 A
Generator Earth Leakage Currents	System meets requirements of IEC 60601-1

Power Cable Specifications

100-120 VAC Installations	
Maximum Cable Length	16.4 ft (5m)
Minimum Rated Cable and Connector Current	10 A
Minimum Rated Cable and Connector Voltage	125 VAC
Cable Connection at Equipment	IEC 60320 type C13 connector

220-240 VAC Installations	
Maximum Cable Length	16.4 ft (5m)
Minimum Rated Cable and Connector Current	5 A
Minimum Rated Cable and Connector Voltage	250 VAC
Cable Connection at Equipment	IEC 60320 type C13 connector

Protection against Electric Shock

Class 1 Equipment

Type CF Applied Parts (device port)

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The equipotential stud on the back panel of the generator is provided for use with other medical electrical systems that require such connections. The equipotential stud is not intended for protective earthing.

RF Energy Output Data

Note: The following load curves apply to the fundamental power delivery capability of the generator alone. They do not imply a given power output for any given device when used with the generator.

Figure 21: RF Energy Output versus Load Impedance – Full Power


Figure 22: RF Energy Output versus Load Impedance – Full Power (Zoomed View)


Figure 23: RF Energy Output versus Load Impedance – Half Power


Figure 24: RF Energy Output versus Load Impedance – Half Power (Zoomed View)


Figure 25: RF Energy Output versus Load Impedance - High Voltage



Figure 26: RF Energy Output versus Load Impedance - High Voltage (Zoomed View)

Electromagnetic Compatibility

As with all medical electrical equipment, the Voyant generator requires special precautions to ensure electromagnetic compatibility (EMC) with other medical electrical devices. To ensure EMC, the Voyant generator must be installed and operated according to the EMC information provided in this manual. The Voyant generator complies with the appropriate IEC 60601-1-2 and 60601-2-2 specifications regarding EMC.

For warnings and precautions related to EMC, see Chapter 2, *Patient & Operating Room Safety*.

Guidance and Manufacturer's Declaration-Electromagnetic Emissions

The Voyant generator is intended for use in the electromagnetic environment specified below. The customer or the user of the Voyant generator should assure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The Voyant generator is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class A	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Guidance and Manufacturer’s Declaration-Electromagnetic Immunity

The Voyant generator is intended for use in the electromagnetic environment specified below. The customer or the user of the Voyant generator should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV Differential Mode ± 2 kV Common Mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions, and voltage variations on power supply input lines IEC 61000-4-11	Voltage dips: 0 % U _r ; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % U _r ; 1 cycle and 70 % U _r ; 25/30 cycles (50/60 Hz) Single phase: at 0° Voltage interruption: 0 % U _r ; 250/300 cycle (50/60 Hz)	Voltage dips: 0 % U _r ; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % U _r ; 1 cycle and 70 % U _r ; 25/30 cycles (50/60 Hz) Single phase: at 0° Voltage interruption: 0 % U _r ; 250/300 cycle (50/60 Hz)	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Voyant generator requires continued operation during power mains interruptions, it is recommended that the Voyant generator be powered from an uninterruptible power supply or a battery.
Power frequency (50-60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note: U_r is the a.c. mains voltage prior to application of the test level.

Guidance and Manufacturer's Declaration-Electromagnetic Immunity

The Voyant generator is intended for use in the electromagnetic environment specified below. The customer or the user of the Voyant generator should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz 6 V in ISM bands between 150 kHz to 80 MHz	3 V	WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the EA020 including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	3 V/m	
Radiated RF Proximity fields from RF wireless communications equipment IEC 61000-4-3	380 MHz to 390 MHz; 27 V/m; PM; 18 Hz 430 MHz to 470 MHz; 28 V/m; (FM ± 5 kHz, 1 kHz sine) 704 MHz to 787 MHz ; 9 V/m; PM; 217 Hz 800 MHz to 960 MHz 28 V/m; PM; 18 Hz 1700 MHz to 1990 MHz 28 V/m; PM; 217 Hz 2400 MHz to 2570 MHz 28 V/m; PM; 217 Hz 5100 MHz to 5800 MHz 9 V/m; PM; 217 Hz		WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the EA020 including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

Note: The emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

Note: The ISM (industrial, scientific and medical) bands between 0.15 MHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz.

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Figura 19:

The screenshot shows a black screen with a green header bar containing the text 'Modo servicio' and a white header bar containing the text 'Resultados 4/10'. Below the headers, there is a table with four columns: 'Rango', 'Resultado', and 'Superado/No superado'. The table contains data for Voltage (V), Current (mA), Power (VA), and Phase. At the bottom of the screen, there are two white buttons: 'Atrás' on the left and 'Continuar' on the right.

	Rango	Resultado	Superado/ No superado
Voltaje (V)	54,81 - 60,58	57,87 57,51	Superado
Corriente (mA)	2879,07 - 3182,13	3030,50 3034,04	Superado
Potencia (VA)	166,12 - 183,60	175,46 174,97	Superado
Fase	(-6,37) - (-2,37)	-5,35 -5,75	Superado

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Figure 19 :

The screenshot shows a black screen with a green border. At the top left, there is a green pill-shaped button labeled 'Mode service'. At the top right, the text 'Résultats 4/10' is displayed. Below this, there is a table with four columns: 'Plage', 'Résultat', and 'Réussite/Échec'. The table contains data for Tension (V), Courant (mA), Puissance (VA), and Phase. At the bottom, there are two white buttons: 'Retour' on the left and 'Continuer' on the right.

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Bild 19:

The screenshot shows a black screen with a green border. At the top left, there is a green pill-shaped button labeled 'Serviceläge'. On the right side, the text 'Resultat 4/10' is displayed. Below this, there is a table with four columns: 'Intervall', 'Resultat', and 'Godkänd/ Icke godkänd'. The first column is empty. At the bottom, there are two small white buttons: 'Tillbaka' on the left and 'Fortsätt' on the right.

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Figur 19:

The screenshot shows a black screen with a green border. At the top left, there is a green pill-shaped button labeled 'Servicemodus'. At the top right, the text 'Resultat 4/10' is displayed. Below this, there is a table with four columns: 'Spenning (V)', 'Strøm (mA)', 'Effekt (VA)', and 'Fase'. Each row contains a range of values for 'Område', a specific value for 'Resultat', and a status for 'Godkjent/feilet'. At the bottom, there are two buttons: 'Tilbake' and 'Fortsett'.

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
	<p>Defibrillation-proof type CF applied part Pieza a prueba de desfibrilación tipo CF Équipement de type CF protégé contre les défibrillateurs Compatibile con defibrillatori con parte applicata di tipo CF Defibrillationssicheres Anwendungssteil des Typs CF Type CF defibrillatiebestending toegepast onderdeel Defibrilleringssikker patientansluten del, typ CF Defibrillationsikker type CF anvendt del Defibrillointisuojattu tyyppi CF liityntäosa Defibrilleringsikker, pasientnær del av type CF Peça aplicada do tipo CF à prova de desfibrilhação</p>	<p>IEC 60601-1:ed3:2005:am1:2012 Medical electrical equipment –Part 1: General requirements for basic safety and essential performance</p>	<p>IEC 60417-5336</p>
	<p>On Activado Marche Acceso Ein Aan På Til Päällä På Ligado</p>	<p>IEC 60601-1:ed3:2005:am1:2012 Medical electrical equipment –Part 1: General requirements for basic safety and essential performance</p>	<p>IEC 60417-5007</p>
	<p>Off Desactivado Arrêt Spento Aus Ult Av Fra Pois päältä Av Desligado</p>	<p>IEC 60601-1:ed3:2005:am1:2012 Medical electrical equipment –Part 1: General requirements for basic safety and essential performance</p>	<p>IEC 60417-5008</p>
	<p>Equipotentiality Equipotencialidad Equipotentialité Equipotenzialità Potenzialgleichheit Equipotential Ekvipotentialitet Ækvipotentialitet Tasapotentialaali Ekvipotentialitet Equipotencialidade</p>	<p>IEC 60601-1:ed3:2005:am1:2012 Medical electrical equipment –Part 1: General requirements for basic safety and essential performance</p>	<p>IEC 60417-5021</p>
	<p>Authorised Representative in the European Community Representante autorizado en la Comunidad Europea Représentant agréé pour la Communauté européenne Rappresentante autorizzato nella Comunità Europea Autorisierte Vertretung in der Europäischen Gemeinschaft Gemechtigde vertegenwoordiger in de Europese Gemeenschap Auktoriserad representant i EU Autoriseret representant i EU Valtuutettu edustaja Euroopan yhteisössä Autorisert representant i Det europeiske fellesskap Representante Autorizado na Comunidade Europeia</p>	<p>ISO 15223-1:2016 Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied — Part 1: General requirements</p>	<p>5.1.2</p>
	<p>Low priority alarm Alarma de prioridad baja Alarme priorité basse Allarme a bassa priorità Alarm geringer Priorität Alarm med lave prioritet Larm av låg prioritet Lavprioritetsalarm Ei tärkeä hälytys Alarm med lav prioritet Alarme de prioridade baixa</p>	<p>IEC 60601-1-8:2006:am1:2012 Medical electrical equipment: Part 1-8: General requirements for basic safety and essential performance – Collateral standard: General requirements, tests, and guidance for alarm systems in medical electrical equipment and medical electrical systems</p>	<p>IEC 60417-5307</p>
	<p>Manufacturer Fabricante Fabricant Produttore Hersteller Fabrikant Tilverkare Producent Valmistaja Produsent Fabricante</p>	<p>ISO 15223-1:2016 Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied — Part 1: General requirements</p>	<p>5.1.1</p>
	<p>Manufacturer and date of manufacture Fabricante y fecha de fabricación Fabricant et date de fabrication Produttore e data di produzione Hersteller und Herstellungsdatum Fabrikant en productiedatum Tilverkare och tillverkningsdatum Producent og produktionsdato Valmistaja ja valmistuspäivämäärä Produsent og produksjonsdato Fabricante e data de fabrico</p>	<p>ISO 15223-1:2016 Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied — Part 1: General requirements</p>	<p>5.1.1</p>

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